

免疫法糞便檢體採集方法

Stool Blood Testing of Immunology Collect Procedure

◎免疫法糞便採集步驟：

1. 檢體清楚的標示姓名與採便的日期、時間
2. 旋開(綠色)瓶蓋，向上拔開取出採便棒。
3. 糞便儘量避免沾到水
4. 如圖示在整條糞便以輕劃的方式在表面不同部位上刮取。刮取量約為完全覆蓋刮取溝槽即可
5. 將綠色蓋子插回插入採檢管後，將蓋子用力壓緊。這個動作只可執行一次，切勿重覆採集。

◎注意事項：

1. 採取大便檢體，採完後放入塑膠袋中，立即送至檢體受理處，否則需放在避光陰涼處。
2. 請勿將管瓶中液體倒出，或任意在瓶中加入水。
3. 遇痔瘡出血及女性月經期間，暫停採便。
4. 採便容器若遺失時，請至檢醫部登記處開立檢體容器醫令後，至掛號處結帳，持收據再回至檢醫部登記處領取採便容器。

*若有相關疑問，請詢問檢驗醫學部(鏡檢暨細胞學組) (02) 23123456 轉分機 267681


◎Collecting methods for stool blood test of Immunology:

1. Please check the accuracy of the specimen label information on the sampling bottle.
2. Screw open the green cap and take out the stirrer.
3. Avoid stool getting into water.
4. Scrape across different parts of the surface of stool. Only a small amount of stool.
5. After inserting the collection tube, firmly press the lid. Once inserted, do not unplug.

◎Announcements:

1. Please collect the stool sample as shown above. After collecting, put it in the plastic bag and take it to SCC in no time. But if you cannot do it right now, please keep it a cool place and out of the light.
2. Please don't pour out the liquid or add water to the sampling bottle.
3. If you have hemorrhoids with bleeding or in menstrual period, please don't collect stool.
4. If you lost the sampling bottle, please go to Laboratory Medicine Registration counter to solve the problem.

*For any questions, please call Laboratory Medicine (Clinical Microscopy and Cytology) for inquiry: (02) 23123456 ext 267681



注意事項 (Precautions) :

- 按照上面圖解採集糞便檢體，完成放入夾鏈袋中並置於陰涼處。
Collect the stool sample as shown above, put it into the zipper bag and keep in a cool place after completed.
 - 勿將管中的液體倒出、或任意於管中加水。
Do not pour out the liquid or add water to the sampling bottle.
 - 遇痔瘡出血或女性月經期間，暫緩採便。
Do not collect sample during bleeding hemorrhoids or menstrual period.
- ※ 請珍惜醫療資源，將採集完畢之容器儘速交回指定醫療院所或衛生所。
Return the collected container to the medical institutions as soon as possible.